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- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3351 - 364 Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref: 1001096	Cont.	R1: 1 x 60 mL, R2: 1 x 20 mL, CAL: 1 x 1 mL
Ref: 1001097		R1: 1 x 30 mL, R2: 1 x 10 mL, CAL: 1 x 1 mL
Ref: 1001098		R1: 1 x 240 mL, R2: 1 x 80 mL, CAL: 1 x 1 mL
Ref: 1001099		R1: 3 x 240 mL, R2: 1 x 240 mL, CAL: 4 x 1 mL

HDL Cholesterol D

Direct. Enzymatic colorimetric

Quantitative determination of HDL cholesterol IVD

Store at 2-8°C

PRINCIPLE OF THE METHOD


Directly determination of serum HDLc (high-density lipoprotein cholesterol) levels without the need for any pre-treatment or centrifugation of the sample^{3,5}.

The assay takes place in two steps.

- 1° Elimination of lipoprotein no-HDL

- 2° Measurement of HDLc

The intensity of the color formed is proportional to the HDLc concentration in the sample.

CLINICAL SIGNIFICANCE

HDL particles are high-density lipoproteins that transport cholesterol from the body tissues to the liver. Since HDL can remove cholesterol from the arteries and carry it back to the liver for their excretion, HDL is known as "good cholesterol" because high levels are thought to lower the risk of heart disease and coronary artery disease.

A low HDL cholesterol levels, is considered a greater heart disease risk^{1,2,4}.

Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

REAGENTS

R 1	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 6,6	100 mM
	N-(2-hydroxy-3-sulfopropyl)-3,5-dimethoxyaniline (HDAOS)	0,7 mM
	Cholesterol Esterase	≥ 800 U/L
	Cholesterol oxidase	≥ 500 U/L
	Catalase	≥ 300 U/L
R 2	N,N-bis(2-hydroxyethyl)-2-aminoethanesulphonic acid pH 7,0	1,1 mmol/L
	4 - Aminoantipyrine (4-AA)	100 mM
	Peroxidase	≥ 3500 U/L
HDLc/ LDLc CAL	Calibrator. Lyophilized human serum.	

PRECAUTIONS

HDLc/ LDLc CAL

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

PREPARATION

- R 1 and R 2: Are ready to use.

- HDLc/ LDLc CAL: Dissolve the contents with 1 mL of distilled water. Cap vial and mix gently to dissolve contents.

STORAGE AND STABILITY ^(Note 1)

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not freeze the reagents.

- HDLc / LDLc CAL: Once reconstitute is stable 30 hours at 20-25°C, 2 weeks at 2-8°C or 3 months at -20°C.

Do not use reagents over the expiration date.

Signs of reagent deterioration:

- Presence of particles and turbidity.

ADDITIONAL EQUIPMENT

- Spectrophotometer or colorimeter measuring at 570 nm.

- Matched cuvettes 1,0 cm light path.

- General laboratory equipment.

SAMPLES

Serum, heparinized plasma or EDTA plasma. If any samples show precipitates, centrifuge before using.

Stability of the sample: 6 days at 2-8°C and 1 year when stored at -70°C.

PROCEDURE

1. Assay conditions:

Wavelength: 550-650 nm

Cuvette: 1 cm light path

Temperature: 37°C

2. Adjust the instrument to zero with distilled water.

3. Pipette into a cuvette:

	Blank	Calibrator	Sample
R 1 (μL)	300	300	300
Calibrator (μL)	--	3	--
Sample (μL)	--	--	3

4. Mix and incubate for 5 min at 37°C and read the absorbance (A1) of the samples and calibrator.

5. Add:

	Blank	Calibrator	Sample
R 2 (μL)	100	100	100

6. Mix and incubate for 5 min at 37°C and read the absorbance (A2) of the samples and calibrator, against the Blank.

CALCULATIONS

$$\frac{(A2-A1) \text{ Sample} - (A2-A1) \text{ Blank}}{(A2-A1) \text{ Calibrator} - (A2-A1) \text{ Blank}} \times \text{Calibrator conc.} = \text{mg/dL of HDL-c in the sample}$$

Conversion factor: mg/dL x 0,0259= mmol/L.

QUALITY CONTROL

Control sera are recommended to monitor the performance of assay procedures: SPINROL H Normal and Pathologic (Ref. 1002120 and 1002210).

If control values are found outside the defined range, check the instrument, reagents and calibrator for problems.

Each laboratory should establish its own Quality Control scheme and corrective actions if controls do not meet the acceptable tolerances.

REFERENCE VALUES

	Men	Women
Low risk	> 50 mg/dL	> 60 mg/dL
Normal risk	35 – 50 mg/dL	45 – 60 mg/dL
High risk	< 35 mg/dL	< 45 mg/dL

These values are for orientation purpose; each laboratory should establish its own reference range.

PERFORMANCE CHARACTERISTICS

Measuring range: From *detection limit* of 5,0 mg/dL to *linearity limit* of 151 mg/dL.

If the results obtained were greater than linearity limit, dilute the sample 1/2 with NaCl 9 g/L and multiply the result by 2.

Precision:

Mean (mg/dL)	Intra-assay (n=20)		Inter-assay (n=20)	
	28,0	76,1	27,5	75,3
SD	0,25	0,81	1,26	2,04
CV (%)	0,89	1,06	4,60	2,71

Sensitivity 1 mg/dL = 0,001399 (A).

Accuracy: Results obtained using SPINREACT reagents (y) did not show systematic differences when compared with other commercial reagents (x).

The results obtained using 50 samples were the following:

Correlation coefficient (r)²: 0,938.

Regression equation: y= 0,9825x - 1,41606.

The results of the performance characteristics depend on the analyzer used.

INTERFERENCES

No interferences were observed to bilirubin up to 30 mg/dL, hemoglobin up to 500 mg/dL, rheumatoid factors up to 1000 IU/mL or lipemia up to 1200 mg/dL.

Lipemic samples with a triglyceride concentration >1200 mg/dL should be diluted 1/10 with NaCl 9 g/L and multiply the result by 10.

NOTES

1. The reagent 2 presents yellowish coloration due to the peroxidase, but it does not affect its functionality.

2. SPINREACT has instruction sheets for several automatic analyzers. Instructions for many of them are available on request.

BIBLIOGRAPHY

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 – 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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BIBLIOGRAPHIE

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3:351 - 364
- Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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BIBLIOGRAFIA

- National Institutes of Health Consensus Development Conference Statement: Triglyceride, High Density Lipoprotein and Coronary Heart Disease. Washington D.C. Feb 26-28, 1992.
- Izawa S., Okada M., Matsui H., and Horita Y. J. Medicine and Pharmaceutical Sci., 1385 - 1388, 37 (1997).
- Shih WJ, Bachorik PS, Haga JA, Myers GL, Stein EA; Clinical Chemistry, 2000; 46:3351 – 364 Third Report of the National Cholesterol Education Programme (NCEP) Expert Panel on Detection, Evaluation and treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA Publication, Vol 285, No. 19, P2486 - 2497; 2001.
- Jacobs, D. et al. In Laboratory and Test Handbook; Jacobs, D.S; Kasten, B.L., De Mott, W.R., Wolfson, W.L., Eds; Lexi - Comp Inc: Hudson (Cleveland), 1990; P. 219.

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Ref: 1001096	Cont.	R1: 1 x 60 mL, R2: 1 x 20 mL, CAL: 1 x 1 mL
Ref: 1001097		R1: 1 x 30 mL, R2: 1 x 10 mL, CAL: 1 x 1 mL
Ref: 1001098		R1: 1 x 240 mL, R2: 1 x 80 mL, CAL: 1 x 1 mL
Ref: 1001099		R1: 3 x 240 mL, R2: 1 x 240 mL, CAL: 4 x 1 mL