SPINCHEM

SPINCHEM 1 / SPINCHEM 2
Human Control serum
Suero Control humano**Serum control multicomponent for clinical chemistry assays IVD**

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

SPINCHEM 1 or SPINCHEM 2 is a human lyophilised serum. With most constituent concentrations and activities in the normal (Spinchem 1) or pathologic (Spinchem 2) range. It's intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Enzymes (human, bovine and porcine). Nonprotein constituents and bacteriostatic agents. The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

Contains human sourced material. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents. Consider all blood derivatives potentially infectious. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical device. Wear suitable gloves and eye/face protection. In case of an accident, seek medical device immediately and show this container or label.

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BIBLIOGRAPHIE

1. Council Directive (2000/54EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.
2. International Federation of Clinical Chemistry (IFCC) Education Division, Expert Panel of Quantities and Units: *A Protocol for the Conversion of Clinical Laboratory data*, Journal of Automatic Chemistry Vol. 11, No 5 (Sept – Oct 1989), pp. 223-226

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SPINCHEM 1 Ref: 1002015 4 x 5 mL
SPINCHEM 2 Ref: 1002025 4 x 5 mL

SPINCHEM

SPINCHEM 1 / SPINCHEM 2
Human Control serum
Soro Controlo humano**Serum control multicomponent for clinical chemistry assays
IVD**

Store at 2 - 8°C.

PRODUCT CHARACTERISTICS

SPINCHEM 1 or SPINCHEM 2 is a human lyophilised serum.
With most constituent concentrations and activities in the normal (Spinchem 1) or pathologic (Spinchem 2) range.
It's intended for control of accuracy for use with manual and automated analytical procedures.

REAGENTS

Human serum. Enzymes (human, bovine and porcine).
Nonprotein constituents and bacteriostatics agents.
The concentration / activities of the components are lot-specific.
The exact values and ranges valid for reagents are given informational purpose in the value sheet.

PRECAUTIONS

Contains human sourced material. No known test method can offer complete assurance that products derived from human blood will not transmit infectious agents. Consider all blood derivatives potentially infectious.
Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical device. Wear suitable gloves and eye/face protection. In case of an accident, seek medical device immediately and show this container or label.

PREPARATION

-Using a Class A volumetric pipette, carefully add the exact amount of distilled or deionized water as specified on the vial label to the control material.
-Swirl the vial gently and allow to stand at room temperature for 5 to 10 minutes.
-Swirl contents until homogeneous. DO NOT VORTEX. Remove the sample required for testing.
Improper handling and/or storage of the control can affect results. Inaccurate reconstitution of the control and errors in assay technique can cause erroneous results.

STORAGE AND STABILITY

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store tightly capped when not use.

- After reconstitution is stable for:
At 2°C to 8°C 10 days
- Direct Bilirubin (stored protected from light):
At 2°C to 8°C 5 days
- Alkaline phosphatase may increase in activity with time.

BIBLIOGRAPHY

1. Council Directive (2000/54EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.
2. International Federation of Clinical Chemistry (IFCC) Education Division, Expert Panel of Quantities and Units: *A Protocol for the Conversion of Clinical Laboratory data*, Journal of Automatic Chemistry Vol. 11, No 5 (Sept – Oct 1989), pp. 223-226

PACKAGING

SPINCHEM 1 Ref: 1002015 Cont. 4 x 5 mL
SPINCHEM 2 Ref: 1002025 Cont. 4 x 5 mL

**Soro controlo multiparamétrico para ensaios de Bioquímica
clínica**IVD
Conservar a 2 - 8°C**CARACTERÍSTICAS DO PRODUTO**

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